



VWA3A rabbit pAb

Cat No.:ES11234

For research use only

Overview

Product Name	VWA3A rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein . at AA range: 710-790
Specificity	VWA3A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	von Willebrand factor A domain-containing protein 3A
Gene Name	VWA3A
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	130kD
Human Gene ID	146177
Human Swiss-Prot Number	A6NCI4
Alternative Names	
Background	similarity:Contains 2 VWFA domains.,

