

## HBA rabbit pAb

## Cat No.:ES11261

For research use only

## Overview

Product Name	HBA rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human;Rat;Mouse	
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000	
Immunogen	Synthesized peptide derived from human protein .	
-	at AA range: 30-110	
Specificity	HBA Polyclonal Antibody detects endogenous levels	
	of protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Hemoglobin subunit alpha (Alpha-globin)	
	(Hemoglobin alpha chain)	
Gene Name	HBA1; HBA2	
Cellular localization	extracellular region, cytosol, hemoglobin	
	complex,membrane,cytosolic small ribosomal	
	subunit, haptoglobin-hemoglobin	
	complex,extracellular exosome,endocytic vesicle	
	lumen,blood microparticle,	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	15kD	
Human Gene ID	3039	
Human Swiss-Prot Number	P69905	
Alternative Names		
Background	The human alpha globin gene cluster located on	
	chromosome 16 spans about 30 kb and includes	
	seven loci: 5'- zeta - pseudozeta - mu -	
	pseudoalpha-1 - alpha-2 - alpha-1 - theta - 3'. The	



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



alpha-2 (HBA2) and alpha-1 (HBA1) coding sequences are identical. These genes differ slightly over the 5' untranslated regions and the introns, but they differ significantly over the 3' untranslated regions. Two alpha chains plus two beta chains constitute HbA, which in normal adult life comprises about 97% of the total hemoglobin; alpha chains combine with delta chains to constitute HbA-2, which with HbF (fetal hemoglobin) makes up the remaining 3% of adult hemoglobin. Alpha thalassemias result from deletions of each of the alpha genes as well as deletions of both HBA2 and HBA1; some nondeletion alpha thalassemias have also been reported. [provided by RefSeq, Jul 2008],



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C