



ELK Biotechnology

SHPS1 Rabbit pAb
Catalog Number: ES11378
For research use only.

Overview

Product Name	SHPS1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	Synthesized peptide derived from human protein . at AA range; 91-140
Specificity	SHPS1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	SHPS1
Gene Name	SIRPA BIT MFR MYD1 PTPNS1 SHPS1 SIRP
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	140885
Human Swiss-Prot Number	P78324
Alternative Names	