



# ELK Biotechnology

**His-Tag Rabbit pAb**  
**Catalog Number: ES1144**  
**For research use only.**

## Overview

---

<b>Product Name</b>	His-Tag polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB</b>
<b>Species Cross-Reactivity</b>	<b>N/A</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/5000</b>
<b>Immunogen</b>	Synthetic Peptide
<b>Specificity</b>	The antibody detects C-terminal, internal, and N-terminal His-tag fusion proteins.
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	His-Tag
<b>Gene Name</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	
<b>Human Swiss-Prot Number</b>	
<b>Alternative Names</b>	