



## ELK Biotechnology

**Plant-actin Rabbit pAb**  
**Catalog Number: ES1152**  
**For research use only.**

### Overview

---

Product Name	Plant-actin polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Arabidopsis Nicotiana tabacum Oryza sativa Hordeum vulgare
Recommended dilutions	Western Blot:1/3000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous plant-actin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Plant-actin
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	43kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	