

## O10P1 rabbit pAb

Cat No.: ES11533

For research use only

## Overview

Product Name O10P1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

**Species Cross-Reactivity** Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

**Specificity** O10P1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Olfactory receptor 10P1 (Olfactory receptor 10P2)

(Olfactory receptor 10P3) (Olfactory receptor

OR12-7)

Gene Name OR10P1 OR10P1P OR10P2P OR10P3P

Cellular localizationCell membrane; Multi-pass membrane protein.PurificationThe antibody was affinity-purified from rabbit<br/>antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 34kD
Human Gene ID 121130
Human Swiss-Prot Number Q8NGE3

**Alternative Names** 

**Background** Olfactory receptors interact with odorant molecules


in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950