



ELK Biotechnology

MeCP2 Rabbit pAb
Catalog Number: ES1155
For research use only.

Overview

Product Name	MeCP2 polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/2000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous MeCP2 proteins.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MeCP2
Gene Name	MECP2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	53kDa
Human Gene ID	4204
Human Swiss-Prot Number	P51608
Alternative Names	