



ELK Biotechnology

Alliinase Rabbit pAb
Catalog Number: ES1160
For research use only.

Overview

Product Name	Alliinase polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Allium sativum L
Recommended dilutions	Western Blot:1/3000
Immunogen	Purified Protein
Specificity	The antibody detects endogenous and recombinant Alliinase proteins.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Alliinase
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	