



ELK Biotechnology

Glucocorticoid receptor Rabbit pAb

Catalog Number: ES1161

For research use only.

Overview

Product Name	Glucocorticoid receptor polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/3000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous GR proteins.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Glucocorticoid receptor
Gene Name	NR3C1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80-97kDa
Human Gene ID	2908
Human Swiss-Prot Number	P04150
Alternative Names	NR3C1; GRL; Glucocorticoid receptor; GR; Nuclear receptor subfamily 3 group C member 1