

## GAPDH Rabbit pAb Catalog Number: ES1162 For research use only.

## Overview

Product Name	GAPDH polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Zebrafish
Recommended dilutions	Western Blot:1/2000-1/5000
Immunogen	Recombinant Protein
Specificity	The antibody detects Zebrafish endogenous GAPDH protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GAPDH
Gene Name	GAPDH
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	37kDa
Human Gene ID	2597
Human Swiss-Prot Number	P04406
Alternative Names	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH