



# ELK Biotechnology

**ECFP Rabbit pAb**

**Catalog Number: ES1163**

**For research use only.**

## Overview

---

Product Name	ECFP polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	N/A
Recommended dilutions	Western Blot:1/2000-1/5000
Immunogen	Recombinant Protein
Specificity	The antibody detects ECFP & CFP tag protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ECFP
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	