



ELK Biotechnology

Lamin B1 Rabbit pAb
Catalog Number: ES1170
For research use only.

Overview

Product Name	Lamin B1 polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/3000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous Lamin B1 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Lamin B1
Gene Name	LMNB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	68,45kDa
Human Gene ID	4001
Human Swiss-Prot Number	P20700
Alternative Names	LMNB1; LMN2; LMNB; Lamin-B1