



ELK Biotechnology

Vinculin Rabbit pAb
Catalog Number: ES1172
For research use only.

Overview

Product Name	Vinculin polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Vinculin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Vinculin
Gene Name	VCL
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	130kDa
Human Gene ID	7414
Human Swiss-Prot Number	P18206
Alternative Names	VCL; Vinculin; Metavinculin