



ELK Biotechnology

H31 Rabbit pAb
Catalog Number: ES11742
For research use only.

Overview

Product Name	H31 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	
Specificity	H31 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	H31
Gene Name	HIST1H3A H3FA; HIST1H3B H3FL; HIST1H3C H3FC; HIST1H3D H3FB; HIST1H3E H3FD; HIST1H3F H3FI; HIST1H3G H3FH; HIST1H3H H3FK; HIST1H3I H3FF; HIST1H3J H3FJ
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14kDa
Human Gene ID	8350
Human Swiss-Prot Number	P68431
Alternative Names	