



# ELK Biotechnology

**HS90A Rabbit pAb**  
**Catalog Number: ES11744**  
**For research use only.**

## Overview

---

Product Name	HS90A polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA</b> <b>1/5000-1/20000</b>
Immunogen	
Specificity	HS90A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HS90A
Gene Name	HSP90AA1 HSP90A HSPC1 HSPCA
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>80kDa</b>
Human Gene ID	3320
Human Swiss-Prot Number	P07900
Alternative Names	