



ELK Biotechnology

α -Skeletal Muscle Actin Rabbit pAb

Catalog Number: ES1176

For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | α -Skeletal Muscle Actin polyclonal antibody |
| Host species | Rabbit |
| Applications | WB |
| Species Cross-Reactivity | Human, Mouse, Rat |
| Recommended dilutions | Western Blot: 1/1000-1/2000 |
| Immunogen | Synthetic Peptide |
| Specificity | The antibody detects endogenous α Skeletal Muscle Actin protein. |
| Formulation | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | α -Skeletal Muscle Actin |
| Gene Name | ACTA1 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 42kDa |
| Human Gene ID | 58 |
| Human Swiss-Prot Number | P68133 |
| Alternative Names | ACTA1; ACTA; Actin, alpha skeletal muscle; Alpha-actin-1 |