



ELK Biotechnology

β II Tubulin Rabbit pAb
Catalog Number: ES1177
For research use only.

Overview

Product Name	β II Tubulin polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human, Mouse, Rat
Recommended dilutions	Western Blot: 1/2000-1/5000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous β II tubulin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	β II Tubulin
Gene Name	TUBB2A
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kDa
Human Gene ID	7280
Human Swiss-Prot Number	Q13885
Alternative Names	