



ELK Biotechnology

Akt3 Rabbit pAb
Catalog Number: ES1181
For research use only.

Overview

Product Name	Akt3 polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous Akt3 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Akt3
Gene Name	AKT3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kDa
Human Gene ID	10000
Human Swiss-Prot Number	Q9Y243
Alternative Names	AKT3; PKBG; RAC-gamma serine/threonine-protein kinase; Protein kinase Akt-3; Protein kinase B gamma; PKB gamma; RAC-PK-gamma; STK-2