



ELK Biotechnology

Pan Methyl Lysine Rabbit pAb
Catalog Number: ES1184
For research use only.

Overview

Product Name	Pan Methyl Lysine polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	All
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous protein Methylated Lysine.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Pan Methyl Lysine
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	