



ELK Biotechnology

IRK6 Rabbit pAb
Catalog Number: ES11847
For research use only.

Overview

Product Name	IRK6 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	
Specificity	IRK6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IRK6
Gene Name	KCNJ6 GIRK2 KATP2 KCNJ7
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	46kDa
Human Gene ID	3763
Human Swiss-Prot Number	P48051
Alternative Names	