



ELK Biotechnology

CK18 Rabbit pAb

Catalog Number: ES1185

For research use only.

Overview

Product Name	CK18 polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous Cytokeratin 18 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CK18
Gene Name	KRT18
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	46kDa
Human Gene ID	3875
Human Swiss-Prot Number	P05783
Alternative Names	KRT18; CYK18; PIG46; Keratin, type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18