



ELK Biotechnology

β Tubulin Rabbit pAb
Catalog Number: ES1186
For research use only.

Overview

Product Name	β Tubulin polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human, Mouse, Rat
Recommended dilutions	Western Blot: 1/2000-1/5000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous β tubulin protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	β Tubulin
Gene Name	TUBB3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	10381
Human Swiss-Prot Number	Q13509
Alternative Names	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III