



ELK Biotechnology

SC6A4 Rabbit pAb
Catalog Number: ES11861
For research use only.

Overview

Product Name	SC6A4 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	
Specificity	SC6A4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	SC6A4
Gene Name	SLC6A4 HTT SERT
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	69kDa
Human Gene ID	6532
Human Swiss-Prot Number	P31645
Alternative Names	