



DLL4 rabbit pAb

Cat No.:ES11956

For research use only

Overview

| | |
|---------------------------------|---|
| Product Name | DLL4 rabbit pAb |
| Host species | Rabbit |
| Applications | WB;ELISA |
| Species Cross-Reactivity | Human;Mouse |
| Recommended dilutions | WB 1:500-2000 ELISA 1:5000-20000 |
| Immunogen | Synthesized peptide derived from part region of human protein |
| Specificity | DLL4 Polyclonal Antibody detects endogenous levels of protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C. Avoid repeated freeze-thaw cycles. |
| Protein Name | Delta-like protein 4 (Drosophila Delta homolog 4) (Delta4) |
| Gene Name | DLL4 UNQ1895/PRO4341 |
| Cellular localization | Cell membrane ; Single-pass type I membrane protein . |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 75kD |
| Human Gene ID | 54567 |
| Human Swiss-Prot Number | Q9NR61 |
| Alternative Names | |
| Background | This gene is a homolog of the Drosophila delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain. [provided by RefSeq, Jul 2008], |

