



## ELK Biotechnology

### Histone H1 (Di Methyl Lys25) Rabbit pAb

Catalog Number: ES1216

For research use only.

#### Overview

---

Product Name	Histone H1 (Di Methyl Lys25) polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human, Mouse, Rat
Recommended dilutions	Western Blot: 1/500-1/1000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Histone H1 (Di Methyl Lys25) protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H1 (Di Methyl Lys25)
Gene Name	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	17-25kDa
Human Gene ID	3024/3006/3007/3008
Human Swiss-Prot Number	Q02539/P16403/P16402/P10412
Alternative Names	HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4