



## ELK Biotechnology

### Histone H2A (Acetyl Lys15) Rabbit pAb

Catalog Number: ES1228

For research use only.

#### Overview

---

Product Name	Histone H2A (Acetyl Lys15) polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human, Mouse, Rat
Recommended dilutions	Western Blot: 1/500-1/1000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Histone H2A (Acetyl Lys15) protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2A (Acetyl Lys15)
Gene Name	HIST1H2AG/HIST1H2AI/HIST1H2AK/HIST1H2AL/HIST1H2AM/HIST2H2A3/HIST2H2AA4/HIST3H2A
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14kDa
Human Gene ID	8329/8330/8332/8336/8969/723790/8337/92815
Human Swiss-Prot Number	P0C0S8/Q6FI13/Q7L7L0
Alternative Names	HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST2H2AA3; H2AFO; HIST2H2AA; HIST2H2AA4; Histone H2A type 2-A; Histone H2A.2; Histone H2A/o; HIST3H2A; Histone H2A type 3