



# ELK Biotechnology

## Histone H2B (Acetyl Lys12) Rabbit pAb

Catalog Number: ES1230

For research use only.

### Overview

---

Product Name	Histone H2B (Acetyl Lys12) polyclonal antibody
Host species	Rabbit
Applications	<b>WB</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/1000-1/2000</b>
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Histone H2B (Acetyl Lys12) protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2B (Acetyl Lys12)
Gene Name	HIST1H2BC
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>14kDa</b>
Human Gene ID	255626/3018/3017/8339/8343/8344/8346/8347
Human Swiss-Prot Number	Q96A08/P33778/P62807
Alternative Names	HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK