



## ELK Biotechnology

### Histone H2A.Z (Acetyl Lys4) Rabbit pAb

Catalog Number: ES1245

For research use only.

#### Overview

---

Product Name	Histone H2A.Z (Acetyl Lys4) polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human, Mouse, Rat
Recommended dilutions	Western Blot: 1/1000-1/2000
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Histone H2A.Z (Acetyl Lys4) protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2A.Z (Acetyl Lys4)
Gene Name	H2AFZ
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14kDa
Human Gene ID	3015
Human Swiss-Prot Number	P0C0S5
Alternative Names	H2AFZ; H2AZ; Histone H2A.Z; H2A/z