

Dok-1 (phospho Tyr362) rabbit pAb

Cat No.:ES1298

For research use only

Overview

Product Name	Dok-1 (phospho Tyr362) rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence:
	1/200 - 1/1000. ELISA: 1/10000. Not yet tested in
	other applications.
Immunogen	The antiserum was produced against synthesized
-	peptide derived from human p62 Dok around the
	phosphorylation site of Tyr362. AA range:329-378
Specificity	Phospho-Dok-1 (Y362) Polyclonal Antibody detects
	endogenous levels of Dok-1 protein only when
	phosphorylated at Y362.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Docking protein 1
Gene Name	DOK1
Cellular localization	[Isoform 1]: Cytoplasm. Nucleus.; [Isoform 3]:
	Cytoplasm, perinuclear region.
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kD
Human Gene ID	1796
Human Swiss-Prot Number	Q99704
Alternative Names	DOK1; Docking protein 1; Downstream of tyrosine
	kinase 1; p62(dok); pp62
Background	docking protein 1(DOK1) Homo sapiens The
	protein encoded by this gene is part of a signal
	transduction pathway downstream of receptor



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C




tyrosine kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multiprotein signaling complexes. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016],

Western Blot analysis of various cells using Phospho-Dok-1 (Y362) Polyclonal Antibody



1.400 1.200 1.000 0.800 0.600 0.400 0.200 0.000 phosphopeptide non-phosphopeptide ■ OD 450nm Reading



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using p62 Dok (Phospho-Tyr362) Antibody

Western blot analysis of lysates from 293 cells, using p62 Dok (Phospho-Tyr362) Antibody. The lane on the right is blocked with the phospho peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C