



ELK Biotechnology

14-3-3 ζ Rabbit pAb
Catalog Number: ES1542
For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | 14-3-3 ζ polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000 |
| Immunogen | Synthesized peptide derived from the Internal region of human 14-3-3 ζ . at AA range; 30-110 |
| Specificity | 14-3-3 ζ Polyclonal Antibody detects endogenous levels of 14-3-3 ζ protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | 14-3-3 ζ |
| Gene Name | YWHAZ |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 30kDa |
| Human Gene ID | 7534 |
| Human Swiss-Prot Number | P63104 |
| Alternative Names | YWHAZ; 14-3-3 protein zeta/delta; Protein kinase C inhibitor protein 1; KCIP-1 |