



ELK Biotechnology

Annexin I Rabbit pAb
Catalog Number: ES1655
For research use only.

Overview

Product Name	Annexin I polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunocytochemistry:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human Annexin I.
Specificity	Annexin I Polyclonal Antibody detects endogenous levels of Annexin I protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Annexin I
Gene Name	ANXA1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	38kDa
Human Gene ID	301
Human Swiss-Prot Number	P04083
Alternative Names	ANXA1; ANX1; LPC1; Annexin A1; Annexin I; Annexin-1; Calpactin II; Calpactin-2; Chromobindin-9; Lipocortin I; Phospholipase A2 inhibitory protein; p35