



ELK Biotechnology

AP-1 Rabbit pAb
Catalog Number: ES1666
For research use only.

Overview

Product Name	AP-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IP,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunoprecipitation:2-5ug/mg lysate ELISA:1/20000
Immunogen	Synthesized peptide derived from human AP-1 around the non-phosphorylation site of Y170.
Specificity	AP-1 Polyclonal Antibody detects endogenous levels of AP-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	AP-1
Gene Name	JUN
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	36kDa
Human Gene ID	3725
Human Swiss-Prot Number	P05412
Alternative Names	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39