



# ELK Biotechnology

**Bcl-2 Rabbit pAb**  
**Catalog Number: ES1759**  
**For research use only.**

## Overview

---

Product Name	Bcl-2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from human Bcl-2 around the non-phosphorylation site of T56.
Specificity	Bcl-2 Polyclonal Antibody detects endogenous levels of Bcl-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bcl-2
Gene Name	BCL2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>26kDa</b>
Human Gene ID	596
Human Swiss-Prot Number	P10415
Alternative Names	BCL2; Apoptosis regulator Bcl-2