



ELK Biotechnology

BRCA1 Rabbit pAb
Catalog Number: ES1787
For research use only.

Overview

Product Name	BRCA1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human BRCA1 around the non-phosphorylation site of S1524.
Specificity	BRCA1 Polyclonal Antibody detects endogenous levels of BRCA1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	BRCA1
Gene Name	BRCA1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	200kDa
Human Gene ID	672
Human Swiss-Prot Number	P38398
Alternative Names	BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53