



# ELK Biotechnology

## Caspase-9 Rabbit pAb Catalog Number: ES1860 For research use only.

### Overview

---

Product Name	Caspase-9 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>Immunofluorescence:1/200-1/1000</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human Caspase-9 around the non-phosphorylation site of S196.
Specificity	Caspase-9 Polyclonal Antibody detects endogenous levels of Caspase-9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Caspase-9
Gene Name	CASP9
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>50kDa</b>
Human Gene ID	842
Human Swiss-Prot Number	P55211
Alternative Names	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6