



## ELK Biotechnology

**Catenin- $\beta$  Rabbit pAb**  
**Catalog Number: ES1863**  
**For research use only.**

### Overview

---

Product Name	Catenin- $\beta$ polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human Catenin- $\beta$ around the non-phosphorylation site of T41/S45.
Specificity	Catenin- $\beta$ Polyclonal Antibody detects endogenous levels of Catenin- $\beta$ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Catenin- $\beta$
Gene Name	CTNNB1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>85kDa</b>
Human Gene ID	1499
Human Swiss-Prot Number	P35222
Alternative Names	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin