



ELK Biotechnology

Caveolin-1 Rabbit pAb Catalog Number: ES1872 For research use only.

Overview

Product Name	Caveolin-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human Caveolin-1 around the non-phosphorylation site of Y14.
Specificity	Caveolin-1 Polyclonal Antibody detects endogenous levels of Caveolin-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Caveolin-1
Gene Name	CAV1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	25kDa
Human Gene ID	857
Human Swiss-Prot Number	Q03135
Alternative Names	CAV1; CAV; Caveolin-1