



ELK Biotechnology

Cbl Rabbit pAb
Catalog Number: ES1877
For research use only.

Overview

Product Name	Cbl polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/5000
Immunogen	Synthesized peptide derived from human Cbl around the non-phosphorylation site of Y700.
Specificity	Cbl Polyclonal Antibody detects endogenous levels of Cbl protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Cbl
Gene Name	CBL
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	100kDa
Human Gene ID	867
Human Swiss-Prot Number	P22681
Alternative Names	CBL; CBL2; RNF55; E3 ubiquitin-protein ligase CBL; Casitas B-lineage lymphoma proto-oncogene; Proto-oncogene c-Cbl; RING finger protein 55; Signal transduction protein CBL