



ELK Biotechnology

CBP Rabbit pAb
Catalog Number: ES1880
For research use only.

Overview

Product Name	CBP polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the C-terminal region of human CBP. at AA range; 2370-2450
Specificity	CBP Polyclonal Antibody detects endogenous levels of CBP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CBP
Gene Name	CREBBP
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	265kDa
Human Gene ID	1387
Human Swiss-Prot Number	Q92793
Alternative Names	CREBBP; CBP; CREB-binding protein