



ELK Biotechnology

Cdk2/Cdc2 Rabbit pAb
Catalog Number: ES1942
For research use only.

Overview

Product Name	Cdk2/Cdc2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from human Cdk2/Cdc2 around the non-phosphorylation site of T160.
Specificity	Cdk2/Cdc2 Polyclonal Antibody detects endogenous levels of Cdk2/Cdc2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Cdk2/Cdc2
Gene Name	CDK2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	34kDa
Human Gene ID	1017
Human Swiss-Prot Number	P24941
Alternative Names	CDK2; CDKN2; Cyclin-dependent kinase 2; Cell division protein kinase 2; p33 protein kinase