



ELK Biotechnology

CYCS Rabbit pAb
Catalog Number: ES2090
For research use only.

Overview

Product Name	CYCS polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat,Monkey
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from the N-terminal region of human CYCS. at AA range; 10-90
Specificity	CYCS Polyclonal Antibody detects endogenous levels of CYCS protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CYCS
Gene Name	CYCS
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15kDa
Human Gene ID	54205
Human Swiss-Prot Number	P99999
Alternative Names	CYCS; CYC; Cytochrome c