



ELK Biotechnology

Dynamin I Rabbit pAb
Catalog Number: ES2208
For research use only.

Overview

Product Name	Dynamin I polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthesized peptide derived from human Dynamin I around the non-phosphorylation site of S774.
Specificity	Dynamin I Polyclonal Antibody detects endogenous levels of Dynamin I protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Dynamin I
Gene Name	DNM1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	97kDa
Human Gene ID	1759
Human Swiss-Prot Number	Q05193
Alternative Names	DNM1; DNM; Dynamin-1