



ELK Biotechnology

Ezrin Rabbit pAb
Catalog Number: ES2311
For research use only.

Overview

Product Name	Ezrin polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human Ezrin around the non-phosphorylation site of Y353.
Specificity	Ezrin Polyclonal Antibody detects endogenous levels of Ezrin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Ezrin
Gene Name	EZR
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	80kDa
Human Gene ID	7430
Human Swiss-Prot Number	P15311
Alternative Names	EZR; VIL2; Ezrin; Cytovillin; Villin-2; p81