

FAM80B rabbit pAb

Cat No.:ES2324

For research use only

Overview

Product Name	FAM80B rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human;Mouse;Rat;Monkey	
Recommended dilutions	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not	
	yet tested in other applications.	
Immunogen	The antiserum was produced against synthesized	
	peptide derived from human RIMKB. AA	
	range:177-226	
Specificity	FAM80B Polyclonal Antibody detects endogenous	
	levels of FAM80B protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.	
Protein Name	Beta-citryl-glutamate synthase B	
Gene Name	RIMKLB	
Cellular localization	Cytoplasm .	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	42kD	
Human Gene ID	57494	
Human Swiss-Prot Number	Q9ULI2	
Alternative Names	RIMKLB; FAM80B; KIAA1238; Beta-citryl-glutamate	
	synthase B; N-acetyl-aspartyl-glutamate synthetase	
	B; NAAG synthetase B; NAAGS; Ribosomal protein S6	
	modification-like protein B	
Background	cofactor:Binds 2 manganese ions per	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,
	subunit., similarity: Belongs to the rimK	
	family., similarity: Contains 1 ATP-grasp domain.,	




+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C







Western blot analysis of lysates from HepG2 and COS cells, using RIMKB Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C