

Fucokinase rabbit pAb

Cat No.:ES2379

For research use only

Overview

Product Name	Fucokinase rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Human;Mouse;Monkey
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence:
	1/200 - 1/1000. ELISA: 1/20000. Not yet tested in
	other applications.
Immunogen	The antiserum was produced against synthesized
	peptide derived from human FUK. AA range:11-60
Specificity	Fucokinase Polyclonal Antibody detects endogenous
	levels of Fucokinase protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	L-fucose kinase
Gene Name	FUK
Cellular localization	cytoplasm, cytosol,
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	110kD
Human Gene ID	197258
Human Swiss-Prot Number	Q8N0W3
Alternative Names	FUK; L-fucose kinase; Fucokinase
Background	fucokinase(FUK) Homo sapiens The protein
	encoded by this gene belongs to the GHMP
	(galacto-, homoserine, mevalonate and
	phosphomevalonate) kinase family and catalyzes the
	phosphorylation of L-fucose to form beta-L-fucose
	1-phosphate. This enzyme catalyzes the first step in
	the utilization of free L-fucose in glycoprotein and



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com


23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



glycolipid synthesis. L-fucose may be important in mediating a number of cell-cell interactions such as blood group antigen recognition, inflammation, and metastatis. While several transcript variants may exist for this gene, the full-length nature of only one has been described to date. [provided by RefSeq, Jul 2008],

Western Blot analysis of various cells using Fucokinase Polyclonal Antibody





COS7

(kD)


170-

130-

95-

72-

55-



Western blot analysis of lysates from COS7 cells, treated with insulin 0.01U/ml 45', using FUK Antibody. The lane on the right is blocked with the synthesized peptide.



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C