



# ELK Biotechnology

**Gab 2 Rabbit pAb**  
**Catalog Number: ES2386**  
**For research use only.**

## Overview

Product Name	Gab 2 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human Gab 2. at AA range; 210-290
Specificity	Gab 2 Polyclonal Antibody detects endogenous levels of Gab 2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Gab 2
Gene Name	GAB2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>75kDa</b>
Human Gene ID	9846
Human Swiss-Prot Number	Q9UQC2
Alternative Names	GAB2; KIAA0571; GRB2-associated-binding protein 2; GRB2-associated binder 2; Growth factor receptor bound protein 2-associated protein 2; pp100