



## ELK Biotechnology

**GABAB R1 Rabbit pAb**  
**Catalog Number: ES2389**  
**For research use only.**

### Overview

---

<b>Product Name</b>	GABAB R1 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,IHC-p,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat,Monkey</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of human GABAB R1. at AA range; 840-920
<b>Specificity</b>	GABAB R1 Polyclonal Antibody detects endogenous levels of GABAB R1 protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	GABAB R1
<b>Gene Name</b>	GABBR1
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>110kDa</b>
<b>Human Gene ID</b>	2550
<b>Human Swiss-Prot Number</b>	Q9UBS5
<b>Alternative Names</b>	GABBR1; GPRC3A; Gamma-aminobutyric acid type B receptor subunit 1; GABA-B receptor 1; GABA-B-R1; GABA-BR1; GABABR1; Gb1