



## ELK Biotechnology

**GATA-1 Rabbit pAb**  
**Catalog Number: ES2402**  
**For research use only.**

### Overview

---

Product Name	GATA-1 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat,Monkey</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/20000</b>
Immunogen	Synthesized peptide derived from human GATA-1 around the non-phosphorylation site of S310.
Specificity	GATA-1 Polyclonal Antibody detects endogenous levels of GATA-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GATA-1
Gene Name	GATA1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>43kDa</b>
Human Gene ID	2623
Human Swiss-Prot Number	P15976
Alternative Names	GATA1; ERYF1; GF1; Erythroid transcription factor; Eryf1; GATA-binding factor 1; GATA-1; GF-1; NF-E1 DNA-binding protein