



ELK Biotechnology

GATA-4 Rabbit pAb
Catalog Number: ES2403
For research use only.

Overview

Product Name	GATA-4 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 Immunofluorescence:1/200-1/1000 ELISA:1/20000
Immunogen	Synthesized peptide derived from human GATA-4 around the non-phosphorylation site of S105.
Specificity	GATA-4 Polyclonal Antibody detects endogenous levels of GATA-4 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GATA-4
Gene Name	GATA4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	48kDa
Human Gene ID	2626
Human Swiss-Prot Number	P43694
Alternative Names	GATA4; Transcription factor GATA-4; GATA-binding factor 4