



ELK Biotechnology

GPR115 Rabbit pAb
Catalog Number: ES2449
For research use only.

Overview

Product Name	GPR115 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/5000
Immunogen	Synthesized peptide derived from the Internal region of human GPR115. at AA range; 300-380
Specificity	GPR115 Polyclonal Antibody detects endogenous levels of GPR115 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GPR115
Gene Name	GPR115
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	78kDa
Human Gene ID	221393
Human Swiss-Prot Number	Q8IZF3
Alternative Names	GPR115; PGR18; Probable G-protein coupled receptor 115; G-protein coupled receptor PGR18